

## OIS WINSTATION 11K™

THE OIS WINSTATION 11K™ DIGITAL IMAGING SYSTEM IS DESIGNED FOR RETINA SPECIALISTS AND RESEARCH AND TEACHING INSTITUTIONS WHO DEMAND THE ABSOLUTE BEST IMAGE QUALITY AVAILABLE. THIS DIGITAL IMAGING SYSTEM USES AN 11 MEGAPIXEL FILM SIZE CCD, RESULTING IN THE MOST DETAILED IMAGES AVAILABLE ON THE MARKET TODAY.

Unsurpassed Image Quality, Ultra-High Resolution





## Highlights

- High-resolution 11 megapixel film size CCD (4008x2672 total pixels)
- Instant FA and color fundus imaging in a single, compact digital camera back
- Seamless integration with Ophthalmology Office™ – EMR
- Reading center export module
- Optional ICG and Autofluorescence

## SYSTEM FEATURES

AutoMontage™  
 AutoImport™  
 Contrast, Sharpen, and Zoom  
 Negative  
 Flip image for laser treatment  
 PDT planning  
 Automatic image overlay  
 Multiple Image Composite  
 Image measurement  
 Macro—batch processing  
 Multiple visit compare  
 Annotation  
 WebStation™ upgrade available for on-line storage and remote access  
 DICOM/VISTA/HL7-compliant

\*Some features are optional upgrades



221 Lathrop Way, Suite I  
 Sacramento, California 95815

Tel: 800.338.8436

Fax: 916.646.0207

www: oisi.com

## Specifications\*\*

Computer	Mid tower Intel® Pentium® processor 1-GB RAM or more Dual 250-GB hard drives or larger 3.5" 1.44MB floppy drive Enhanced USB keyboard Two-button USB scroll mouse UPS battery backup
Monitor	21" Flat panel high-resolution color
Archiving	16X/48X DVD-recorder and CD-recorder
Networking and Telemedicine	100/1000 Mb Ethernet 56K V92 modem Remote access software
Image Capture	11 megapixel (4008x2672 total pixels), high-resolution B&W and color CCD camera, optical interface for fundus camera, digital sync controller and interface cabling. Optional 12-Bit CCD for ICG capabilities.
Software	Microsoft Windows® OIS WinStation™ PC Anywhere® Norton AntiVirus*
Dimensions	21" x 43" footprint 100V/50-60 Hz - 240V/ 50-60 Hz
Warranty	1 year parts and labor (extended warranty available)
Accessories	Equipment table Equipment cart Printer cart Printer Stereo Viewer

\*All brand names are the property of their respective owners.

\*\*OIS strives to continuously update and improve upon our technology. Consequently, product specifications may change without notice. Image resolution will vary depending on fundus camera and video adapter in use.