


OIS WINSTATION 5000™

THE OIS WINSTATION 5000™ DIGITAL IMAGING SYSTEM IS DESIGNED TO MEET THE MORE DEMANDING NEEDS OF RETINA SPECIALISTS WHO REQUIRE THE ABSOLUTE BEST IMAGE QUALITY. UNMATCHED IN THE EYECARE INDUSTRY, THE OIS WINSTATION 5000™ USES A 6 MEGAPIXEL FILM SIZE CCD, RESULTING IN IMAGES THAT SURPASS THE QUALITY OF PHOTOGRAPHIC FILM. THE OIS WINSTATION 5000 IS APPROVED FOR CLINICAL RESEARCH STUDIES.

Crystal-Clear Resolution





SYSTEM FEATURES

AutoMontage™
 AutoImport™
 Contrast, Sharpen, and Zoom
 Negative
 Flip image for laser treatment
 PDT planning
 Automatic image overlay
 Multiple Image Composite
 Image measurement
 Macro—batch processing
 Multiple visit compare
 Annotation
 DICOM-compliant

*Some features are optional upgrades



221 Lathrop Way, Suite I
 Sacramento, California 95815

Tel: 800.338.8436

Fax: 916.646.0207

www: oisi.com

Highlights

- Approved for clinical studies
- Ultra-high 6 megapixel film size CCD (2400x2048 actual image size) resolution that rivals photographic film
- Instant FA and color fundus imaging in a single, compact digital camera back
- Film port mountable to most popular fundus cameras, including Zeiss FF3/4
- Reading center export module
- Optional ICG and Autofluorescence

Specifications**

Computer	Mid tower Intel® Pentium® processor 1-GB RAM or more Dual 250-GB hard drives or larger 3.5" 1.44-MB floppy drive Enhanced keyboard Two-Button Scroll mouse UPS battery backup
Monitor	21" flat panel color monitor
Archiving	16X/48X DVD-recorder and CD-recorder
Networking and Telemedicine	100/1000 Mb Ethernet 56K V92 modem Remote access software
Image Capture	6 megapixel (2400x2048 actual image size), ultra high-resolution B&W and color film-size CCD camera, optical interface for fundus camera, digital sync controller, and interface cabling. Optional 12-Bit CCD for ICG capabilities.
Software	Microsoft Windows® OIS WinStation™ PC Anywhere® Norton AntiVirus®
Dimensions	21"x43" footprint 100V/50-60 Hz–240V/50-60 Hz
Warranty	1 year parts and labor (extended warranty available)
Accessories	Equipment table Equipment cart Printer cart Printer Stereo viewer

*All brand names are the property of their respective owners.

**OIS strives to continuously update and improve upon our technology. Consequently, product specifications may change without notice. Image resolution will vary depending on fundus camera and video adapter in use.